




Product specification

Pipette tip, 1,250 µl, transparent, 384 piece(s)/bag



Product information

Product number	70.1186
----------------	---------

Product descriptions

Name	Pipette tip, 1,250 µl, transparent, 384 piece(s)/bag
Product number	70.1186
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 384 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Length of product	103 mm
Max. volume	1250 µl

Product characteristics

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order quantity	3840
Piece(s) / inner box	384
Piece(s) / case	3840
Piece(s) / pallet	134400
Height of case	443 mm
Width of case	258 mm
Depth of case	347 mm
Volume of case	0.0397 cbm
Weight of case	4.5696 kg
EAN inner box	4038917225236
EAN case	4038917225229

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT