




## Product specification

**Pipette tip, 1,250 µl, transparent, 384 piece(s)/bag**



### Product information

Product number	70.1186
----------------	---------

### Product descriptions

Name	Pipette tip, 1,250 µl, transparent, 384 piece(s)/bag
Product number	70.1186
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 384 piece(s)/bag

### Material & colours

Color of product	transparent
------------------	-------------

### Size

Length of product	103 mm
Max. volume	1250 µl

### Product characteristics

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

### Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Piece(s) / inner box	384
Piece(s) / case	3840
Piece(s) / pallet	134400
Height of inner box	40 mm
Width of inner box	280 mm
Depth of inner box	610 mm
Height of case	260 mm
Width of case	350 mm
Depth of case	445 mm
Volume of case	0.0405 cbm
Weight of case	4.5696 kg
EAN inner box	4038917225236
EAN case	4038917225229

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT