




Product specification

Pipette tip, 1,250 µl, transparent, Biosphere®, 96 piece(s)/box



Product information

Product number	70.1186.200
----------------	-------------

Product descriptions

Name	Pipette tip, 1,250 µl, transparent, Biosphere®, 96 piece(s)/box
Product number	70.1186.200
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, Biosphere®, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 96 piece(s)/box

Material & colours

Color of product	transparent
------------------	-------------

Size

Length of product	103 mm
Max. volume	1250 µl

Product characteristics

Filter	no
Graduation	yes

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	Biosphere
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	96
Piece(s) / case	768
Piece(s) / pallet	79872
Height of inner box	115 mm
Width of inner box	95 mm
Depth of inner box	125 mm
Height of case	121 mm
Width of case	270 mm
Depth of case	413 mm
Volume of case	0.0135 cbm
Weight of case	1.9999 kg
EAN inner box	4038917160377
EAN case	4038917160360

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT