




# Product specification

**Filter tip, 1,250 µl, transparent, Biosphere®, 96 piece(s)/box**



## Product information

Product number	70.1186.210
----------------	-------------

## Product descriptions

Name	Filter tip, 1,250 µl, transparent, Biosphere®, 96 piece(s)/box
Product number	70.1186.210
Product description	Filter tip, working volume: 1,250 µl, transparent, filling level rings, Biosphere®, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 96 piece(s)/box

## Material & colours

Color of product	transparent
Filter material	Polyethylene (PE)

## Size

Length of product	103 mm
Max. volume	1250 µl

## Product characteristics

Filter	yes
Graduation	yes

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	Biosphere
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Piece(s) / inner box	96
Piece(s) / case	768
Piece(s) / pallet	90000
Height of inner box	115 mm
Width of inner box	95 mm
Depth of inner box	125 mm
Height of case	121 mm
Width of case	270 mm
Depth of case	413 mm
Volume of case	0.0135 cbm
Weight of case	2.0797 kg
EAN inner box	4038917126496
EAN case	4038917080354

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT