




Product specification

Pipette tip, 10 ml, transparent, 100 piece(s)/bag



Product information

Product number	70.1187.001
----------------	-------------

Product descriptions

Name	Pipette tip, 10 ml, transparent, 100 piece(s)/bag
Product number	70.1187.001
Product description	Pipette tip, working volume: 10 ml, transparent, filling level rings, suitable for Eppendorf, Gilson and Brand and products of identical design, 100 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Length of product	160.8 mm
Max. volume	10 ml

Product characteristics

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	20000
Height of inner box	40 mm
Width of inner box	325 mm
Depth of inner box	520 mm
Height of case	400 mm
Width of case	289 mm
Depth of case	467 mm
Volume of case	0.054 cbm
Weight of case	6.6 kg
EAN inner box	4038917574617
EAN case	4038917574600

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT