




Product specification

Pipette tip, extra long, 200 µl, transparent, 96 piece(s)/box



Product information

Product number	70.1189.105
----------------	-------------

Product descriptions

Name	Pipette tip, extra long, 200 µl, transparent, 96 piece(s)/box
Product number	70.1189.105
Product description	Pipette tip, working volume: 200 µl, transparent, filling level rings, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and products of identical design, extra long, 96 piece(s)/box

Material & colours

Color of product	transparent
------------------	-------------

Size

Length of product	77 mm
Max. volume	200 µl

Product characteristics

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Piece(s) / inner box	96
Piece(s) / case	960
Piece(s) / pallet	120960
Height of inner box	85 mm
Width of inner box	84 mm
Depth of inner box	120 mm
Height of case	105 mm
Width of case	250 mm
Depth of case	420 mm
Volume of case	0.011 cbm
Weight of case	1.9402 kg
EAN inner box	4038917314497
EAN case	4038917314480

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT