



Product specification

Pipette tip, 1,000 µl, black, 192 piece(s)/folding box



Product information

Product number	70.1191.103
----------------	-------------

Product descriptions

Name	Pipette tip, 1,000 µl, black, 192 piece(s)/folding box
Product number	70.1191.103
Product description	Pipette tip, working volume: 1,000 µl, black, suitable for Application on Sarstedt aliquoters, 192 piece(s)/folding box

Material & colours

Color of product	black
------------------	-------

Size

Max. volume	1000 µl
-------------	---------

Product characteristics

Filter	yes
--------	-----

Purity & certification

Product category	no medical device IVD
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	960
Piece(s) / case	9600
Piece(s) / pallet	115200
Height of case	300 mm
Width of case	480 mm
Depth of case	720 mm
Volume of case	0.1037 cbm
Weight of case	16.3488 kg
EAN inner box	4038917359115
EAN case	4038917381123

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT