


Product specification

Filter tip, 20 µl, transparent, Biosphere® plus, 96 piece(s)/box



Product information

Product number	70.760.213
----------------	------------

Product descriptions

Name	Filter tip, 20 µl, transparent, Biosphere® plus, 96 piece(s)/Box
------	--

Product number	70.760.213
----------------	------------

Product description	Filter tip, working volume: 20 µl, transparent, filling level rings, Biosphere® plus, suitable for Eppendorf, Gilson and Finnpiquette and products of identical design, 96 piece(s)/box
---------------------	---

Product characteristics

Filter	yes
--------	-----

Graduation	yes
------------	-----

Size

Length of product	51 mm
-------------------	-------

Max. volume	20 µl
-------------	-------

Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

Color of product	transparent
------------------	-------------

Filter material	Polyethylene (PE)
-----------------	-------------------

Purity & certification

Product category	no medical device IVD
------------------	-------------------------

Purity standard	Biosphere plus
-----------------	----------------

Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	157440
Height of inner box	131 mm
Width of inner box	89 mm
Depth of inner box	314 mm
Height of case	150 mm
Width of case	326 mm
Depth of case	370 mm
Volume of case	0.0181 cbm
Weight of case	3.2006 kg
EAN inner box	4038917120432
EAN case	4038917017282

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

