



# Product specification

**PCR film, free from DNase/RNase, material: PO,  
transparent**



## Product information

Product number	95.1994
----------------	---------

## Product descriptions

Name	PCR film, free from DNase/RNase, material: PO, transparent
Product number	95.1994
Product description	PCR film, free from DNase/RNase, material: PO, transparent, suitable for qPCR, self-adhesive, application area: -70 - 110 °C, strong adhesion, high chemical resistance, 100 piece(s)/bag

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Width of product	79 mm
Length of product	141 mm

## Product characteristics

Light protection	no
Closure type	Film
Pierceable	no

Suitable for qPCR	yes
Minimum temperature	-70 °C
Maximum temperature	110 °C

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	free from DNase/RNase
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1200
Piece(s) / pallet	288000
Height of inner box	27 mm
Width of inner box	87 mm
Depth of inner box	150 mm
Height of case	160 mm
Width of case	180 mm
Depth of case	180 mm
Volume of case	0.0052 cbm
Weight of case	2.7204 kg
EAN inner box	4038917088114
EAN case	4038917044301

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT