




## Product specification

**PCR plate half skirt, 96 well, transparent, High Profile, 200 µl, PCR Performance Tested, PP**



## Product information

Product number	72.1979.102
----------------	-------------

## Product descriptions

Name	PCR plate half skirt, 96 well, transparent, High Profile, 200 µl, PCR Performance Tested, PP
Product number	72.1979.102
Product description	PCR plate half skirt, 96 well, transparent, High Profile, 200 µl, PCR Performance Tested, material: PP, flat deck, 5 piece(s)/minigrip bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Profile	High Profile

## Size

Diameter	4.95 mm
Width of product	83.6 mm
Length of product	123.6 mm
Height of product	19.33 mm
Max. volume	200 µl

Sample volume

0.2 ml

## Product characteristics

Number of wells	96
Construction of the edges of the plate	half skirt
Base shape	conical
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	3000 U/min
Centrifugation temperature	22 °C
Centrifugation time	5 Min

## Purity & certification

Product category	no medical device   IVD
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Minimum order quantity	100
Piece(s) / inner box	25
Piece(s) / case	100
Piece(s) / pallet	6000
Height of case	175 mm
Width of case	358 mm
Depth of case	398 mm
Volume of case	0.0249 cbm
Weight of case	2.63 kg
EAN inner box	4038917307253
EAN case	4038917307246

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



# SARSTEDT

This document was prepared by EDP support and is valid without signature.