



Product specification

**PCR film, free from DNase/RNase, material: PL,
transparent**



Product information

Product number	95.1993
----------------	---------

Product descriptions

Name	PCR film, free from DNase/RNase, material: PL, transparent
Product number	95.1993
Product description	PCR film, free from DNase/RNase, material: PL, transparent, suitable for qPCR, self-adhesive, application area: -40 - 120 °C, thin material, high optical clarity, 100 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Width of product	80 mm
Length of product	137 mm

Product characteristics

Light protection	no
Closure type	Film
Pierceable	no

Suitable for qPCR	yes
Minimum temperature	-40 °C
Maximum temperature	120 °C

Purity & certification

Product category	no medical device IVD
Purity standard	free from DNase/RNase
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1600
Piece(s) / pallet	192000
Height of inner box	116 mm
Width of inner box	148 mm
Depth of inner box	148 mm
Height of case	203 mm
Width of case	216 mm
Depth of case	216 mm
Volume of case	0.0095 cbm
Weight of case	3.7408 kg
EAN inner box	4038917303569
EAN case	4038917388221

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT