



Product specification

CryoPure tubes, 1.2 ml, QuickSeal screw cap, green



Product information

| | |
|----------------|------------|
| Product number | 72.377.005 |
|----------------|------------|

Product descriptions

| | |
|---------------------|---|
| Name | CryoPure tubes, 1.2 ml, Quickseal screw cap, green |
| Product number | 72.377.005 |
| Product description | CryoPure tubes, 1.2 ml, tube: PP, QuickSeal screw cap, cap assembled, HD-PE, green, external thread, nominal volume: 1 ml, Cryo Performance Tested, 50 piece(s)/bag |

Material & colours

| | |
|------------------|-----------------------------------|
| Tube material | Polypropylene (PP) |
| Closure material | High Density Polyethylene (HD-PE) |
| Color of cap | green |

Size

| | |
|----------------------|----------|
| Diameter | 12.65 mm |
| Length excluding cap | 29.8 mm |
| Length of product | 36.7 mm |
| Max. volume | 1.2 ml |
| Sample volume | 1 ml |

Product characteristics

| | |
|----------------------|---------------------|
| Label/Print | With print |
| Color of print/label | white |
| Vessel type | Tube |
| Thread | External thread |
| Graduation | yes |
| Cap | Quickseal screw cap |
| Closure type | Screw cap |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | diagnostic device, CE |
| CE certified | CE |
| Purity standard | Cryo Performance Tested |
| Sterilisation | Electron irradiation |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|------------------------|---------------|
| Minimum order quantity | 500 |
| Piece(s) / inner box | 500 |
| Piece(s) / case | 2000 |
| Piece(s) / pallet | 104000 |
| Height of inner box | 118 mm |
| Width of inner box | 183 mm |
| Depth of inner box | 293 mm |
| Height of case | 135 mm |
| Width of case | 378 mm |
| Depth of case | 598 mm |
| Volume of case | 0.0305 cbm |
| Weight of case | 4.26 kg |
| EAN inner box | 4038917106016 |
| EAN case | 4038917080538 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT