



Product specification

CryoPure tubes, 2 ml, QuickSeal screw cap, white



Product information

| | |
|----------------|--------|
| Product number | 72.380 |
|----------------|--------|

Product descriptions

| | |
|---------------------|---|
| Name | CryoPure tubes, 2 ml, Quickseal screw cap, white |
| Product number | 72.380 |
| Product description | CryoPure tubes, 2 ml, tube: PP, QuickSeal screw cap, cap assembled, PP, white, Internal thread, with silicon O-ring, nominal volume: 1.6 ml, Cryo Performance Tested, 50 piece(s)/bag |

Material & colours

| | |
|------------------|--------------------|
| Tube material | Polypropylene (PP) |
| O-Ring material | silicone |
| Closure material | Polypropylene (PP) |
| Color of cap | white |

Size

| | |
|----------------------|----------|
| Diameter | 12.65 mm |
| Length excluding cap | 37.8 mm |
| Length of product | 47.7 mm |
| Max. volume | 2 ml |
| Sample volume | 1.6 ml |

Product characteristics

| | |
|----------------------|---------------------|
| Design | with silicon O-ring |
| Label/Print | With print |
| Color of print/label | white |
| Vessel type | Tube |
| Thread | Internal thread |
| Graduation | yes |
| Cap | Quickseal screw cap |
| Closure type | Screw cap |

Purity & certification

| | |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Purity standard | Cryo Performance Tested |
| Sterilisation | Electron irradiation |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 500 |
| Piece(s) / case | 2000 |
| Piece(s) / pallet | 104000 |
| Height of inner box | 118 mm |
| Width of inner box | 183 mm |
| Depth of inner box | 293 mm |
| Height of case | 135 mm |
| Width of case | 378 mm |
| Depth of case | 598 mm |
| Volume of case | 0.0305 cbm |
| Weight of case | 4.93 kg |
| EAN inner box | 4038917081405 |
| EAN case | 4038917080606 |

This document was prepared by EDP support
and is valid without signature.



SARSTEDT