



Product specification

CryoPure tubes, 2 ml, QuickSeal screw cap, red



Product information

Product number	72.380.002
----------------	------------

Product descriptions

Name	CryoPure tubes, 2 ml, Quickseal screw cap, red
Product number	72.380.002
Product description	CryoPure tubes, 2 ml, tube: PP, QuickSeal screw cap, cap assembled, PP, red, Internal thread, with silicon O-ring, nominal volume: 1.6 ml, Cryo Performance Tested, 50 piece(s)/bag

Material & colours

Tube material	Polypropylene (PP)
O-Ring material	silicone
Closure material	Polypropylene (PP)
Color of cap	red

Size

Diameter	12.65 mm
Length excluding cap	37.8 mm
Length of product	47.7 mm
Max. volume	2 ml
Sample volume	1.6 ml

Product characteristics

Design	with silicon O-ring
Label/Print	With print
Color of print/label	white
Vessel type	Tube
Thread	Internal thread
Graduation	yes
Cap	Quickseal screw cap
Closure type	Screw cap

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	Cryo Performance Tested
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	104000
Height of inner box	118 mm
Width of inner box	183 mm
Depth of inner box	293 mm
Height of case	135 mm
Width of case	378 mm
Depth of case	598 mm
Volume of case	0.0305 cbm
Weight of case	4.97 kg
EAN inner box	4038917106078
EAN case	4038917080613

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT