



Product specification

CryoPure tubes, 2 ml, QuickSeal screw cap, yellow



Product information

Product number	72.380.004
----------------	------------

Product descriptions

Name	CryoPure tubes, 2 ml, Quickseal screw cap, yellow
Product number	72.380.004
Product description	CryoPure tubes, 2 ml, tube: PP, QuickSeal screw cap, cap assembled, PP, yellow, Internal thread, with silicon O-ring, nominal volume: 1.6 ml, Cryo Performance Tested, 50 piece(s)/bag

Material & colours

Tube material	Polypropylene (PP)
O-Ring material	silicone
Closure material	Polypropylene (PP)
Color of cap	yellow

Size

Diameter	12.65 mm
Length excluding cap	37.8 mm
Length of product	47.7 mm
Max. volume	2 ml
Sample volume	1.6 ml

Product characteristics

Design	with silicon O-ring
Label/Print	With print
Color of print/label	white
Vessel type	Tube
Thread	Internal thread
Graduation	yes
Cap	Quickseal screw cap
Closure type	Screw cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Cryo Performance Tested
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	104000
Height of inner box	118 mm
Width of inner box	183 mm
Depth of inner box	293 mm
Height of case	135 mm
Width of case	378 mm
Depth of case	598 mm
Volume of case	0.0305 cbm
Weight of case	4.96 kg
EAN inner box	4038917106085
EAN case	4038917080620

This document was prepared by EDP support
and is valid without signature.



SARSTEDT