

Product specification

CryoPure tubes, 2 ml, QuickSeal screw cap, yellow



Product information

Product number	72.380.004
----------------	------------

Product descriptions

Name	CryoPure tubes, 2 ml, Quickseal screw cap, yellow
Product number	72.380.004
Product description	CryoPure tubes, 2 ml, tube: PP, QuickSeal screw cap, HD-PE, yellow, Internal thread, with silicon O-ring, nominal volume: 1.6 ml, Cryo Performance Tested, 50 piece(s)/bag

Product characteristics

Design	with silicon O-ring
Label/Print	With print
Color of print/label	white
Tube	Tube
Thread	Internal thread
Graduation	yes
Cap	Quickseal screw cap
Type of cap	Screw cap

Size

Diameter	12.65 mm
Length excluding cap	37.8 mm
Length of product	47.7 mm
Max. volume	2 ml
Sample volume	1.6 ml

Material & colours

Tube material	Polypropylene (PP)
---------------	--------------------

O-Ring material	silicone
Cap material	High Density Polyethylene (HD-PE)
Color of cap	yellow

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	Cryo Performance Tested
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	104000
Height of inner box	118 mm
Width of inner box	183 mm
Depth of inner box	293 mm
Height of case	135 mm
Width of case	378 mm
Depth of case	598 mm
Volume of case	0.0305 cbm
Weight of case	4.96 kg
EAN inner box	4038917106085
EAN case	4038917080620

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

