



## Product specification

### CryoPure tubes, 5 ml, QuickSeal screw cap, red



## Product information

Product number	72.383.002
----------------	------------

## Product descriptions

Name	CryoPure tubes, 5 ml, Quickseal screw cap, red
Product number	72.383.002
Product description	CryoPure tubes, 5 ml, tube: PP, QuickSeal screw cap, cap assembled, HD-PE, red, external thread, nominal volume: 4.5 ml, Cryo Performance Tested, 25 piece(s)/bag

## Material & colours

Tube material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	red

## Size

Diameter	12.65 mm
Length excluding cap	83.8 mm
Length of product	90.7 mm
Max. volume	5 ml
Sample volume	4.5 ml

## Product characteristics

Label/Print	With print
Color of print/label	white
Vessel type	Tube
Thread	External thread
Graduation	yes
Cap	Quickseal screw cap
Closure type	Screw cap

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	Cryo Performance Tested
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

## Packaging

Minimum order quantity	250
Piece(s) / inner box	250
Piece(s) / case	1000
Piece(s) / pallet	52000
Height of inner box	118 mm
Width of inner box	183 mm
Depth of inner box	293 mm
Height of case	135 mm
Width of case	378 mm
Depth of case	598 mm
Volume of case	0.0305 cbm
Weight of case	5.05 kg
EAN inner box	4038917106115
EAN case	4038917080668

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support  
and is valid without signature.



# SARSTEDT