



## Product specification

### Screw cap micro tube, 1.5 ml



## Product information

Product number	72.607
----------------	--------

## Product descriptions

Name	Screw cap micro tube, 1.5 ml
Product number	72.607
Product description	Screw cap micro tube, working volume: 1.5 ml, conical base, with knurling, transparent, without cap, 500 piece(s)/bag

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Diameter	10.8 mm
Length of product	44 mm

## Product characteristics

Base shape	Conical base
Graduation	no
Knurls	yes
Centrifugation max (RCF)	20000 x g

## Purity & certification

Satisfies the requirement	IATA
---------------------------	------

Product category

In-vitro diagnostic (IVD)

CE certified

CE - manufacturer's self-declaration

Batched

yes

## Packaging

Piece(s) / inner box

500

Piece(s) / case

5000

Piece(s) / pallet

210000

Height of inner box

50 mm

Width of inner box

315 mm

Depth of inner box

335 mm

Height of case

250 mm

Width of case

378 mm

Depth of case

375 mm

Volume of case

0.0354 cbm

Weight of case

5.87 kg

EAN inner box

4038917082839

EAN case

4038917001854

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT