



Product specification

Screw cap micro tube, 2 ml



Product information

Product number	72.609.300
----------------	------------

Product descriptions

Name	Screw cap micro tube, 2 ml
Product number	72.609.300
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, without cap, 500 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Diameter	10.8 mm
Length of product	44 mm

Product characteristics

Base shape	skirted conical base
Graduation	no
Knurls	yes
Centrifugation max (RCF)	20000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	180000
Volume of case	0.0492 cbm
Weight of case	6.83 kg
EAN inner box	4038917434287
EAN case	4038917434270

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT