



## Product specification

### Screw cap micro tube, 2 ml



## Product information

Product number	72.609.300
----------------	------------

## Product descriptions

Name	Screw cap micro tube, 2 ml
Product number	72.609.300
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, without cap, 500 piece(s)/bag

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Diameter	10.8 mm
Length of product	44 mm

## Product characteristics

Base shape	skirted conical base
Graduation	no
Knurls	yes
Centrifugation max (RCF)	20000 x g

## Purity & certification

Satisfies the requirement	IATA
---------------------------	------

Product category	In-vitro diagnostic (IVD)
------------------	---------------------------

CE certified	CE
--------------	----

Batched	yes
---------	-----

## Packaging

Piece(s) / inner box	500
----------------------	-----

Piece(s) / case	5000
-----------------	------

Piece(s) / pallet	180000
-------------------	--------

Height of case	246 mm
----------------	--------

Width of case	450 mm
---------------	--------

Depth of case	450 mm
---------------	--------

Volume of case	0.0498 cbm
----------------	------------

Weight of case	7.056 kg
----------------	----------

EAN inner box	4038917434287
---------------	---------------

EAN case	4038917434270
----------	---------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT