



Product specification

Screw cap micro tube, 2 ml, sterile



Product information

Product number	72.609.311
----------------	------------

Product descriptions

Name	Screw cap micro tube, 2 ml, sterile
Product number	72.609.311
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, without cap, sterile, 500 piece(s)/double bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	10.8 mm
Length of product	44 mm

Product characteristics

Base shape	skirted conical base
Graduation	no
Knurls	yes
Centrifugation max (RCF)	20000 x g

Purity & certification

Satisfies the requirement	IATA
---------------------------	------

Product category	In-vitro diagnostic (IVD)
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	2500
Piece(s) / pallet	140000
Volume of case	0.029 cbm
Weight of case	3.63 kg
EAN inner box	4038917434379
EAN case	4038917434362

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT