

Product specification

Screw cap micro tube, 1.5 ml



Product information

Product number	72.687
----------------	--------

Product descriptions

Name	Screw cap micro tube, 1.5 ml
Product number	72.687
Product description	Screw cap micro tube, working volume: 1.5 ml, conical base, transparent, without cap, 500 piece(s)/bag

Product characteristics

Base shape	Conical base
Graduation	no
Knurls	no
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	44 mm

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	yes

Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	210000
Height of case	250 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0367 cbm
Weight of case	5.92 kg
EAN inner box	4038917083393
EAN case	4038917002318

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

