



Product specification

Reaction tube, 1.5 ml, PP



Product information

Product number	72.688
----------------	--------

Product descriptions

Name	Reaction tube, 1.5 ml, PP
Product number	72.688
Product description	Reaction tube, working volume: 1.5 ml, material: PP, transparent, cap: transparent, retaining loop cap, cap attached, 500 piece(s)/bag

Material & colours

Color of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	transparent

Size

Diameter	10.8 mm
Length of product	39 mm

Product characteristics

Graduation	no
Cap	Retaining loop cap
Centrifugation max (RCF)	20000 x g

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	180000
Height of case	310 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0491 cbm
Weight of case	6.59 kg
EAN inner box	4038917083409
EAN case	4038917002325

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT