



## Product specification

### Reaction tube, 2 ml, PP



## Product information

|                |            |
|----------------|------------|
| Product number | 72.689.003 |
|----------------|------------|

## Product descriptions

|                     |  |
|---------------------|--|
| Name                | Reaction tube, 2 ml, PP  |
| Product number      | 72.689.003   |
| Product description | Reaction tube, working volume: 2 ml, material: PP, transparent, cap: red, retaining loop cap, cap attached, 500 piece(s)/bag |

## Material & colours

|                  |                                  |
|------------------|----------------------------------|
| Color of product | transparent                      |
| Closure material | Low Density Polyethylene (LD-PE) |
| Color of cap     | red                              |

## Size

|                   |         |
|-------------------|---------|
| Diameter          | 10.8 mm |
| Length of product | 39 mm   |

## Product characteristics

|                          |                    |
|--------------------------|--------------------|
| Graduation               | no                 |
| Cap                      | Retaining loop cap |
| Centrifugation max (RCF) | 20000 x g          |

## Purity & certification

|                  |                                      |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD)            |
| CE certified     | CE - manufacturer's self-declaration |
| Batched          | yes                                  |

## Packaging

|                      |               |
|----------------------|---------------|
| Piece(s) / inner box | 500           |
| Piece(s) / case      | 5000          |
| Piece(s) / pallet    | 150000        |
| Height of case       | 340 mm        |
| Width of case        | 398 mm        |
| Depth of case        | 398 mm        |
| Volume of case       | 0.0539 cbm    |
| Weight of case       | 7.48 kg       |
| EAN inner box        | 4038917116442 |
| EAN case             | 4038917015134 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT