



Product specification

Reaction tube, 2 ml, PP



Product information

Product number	72.689.004
----------------	------------

Product descriptions

Name	Reaction tube, 2 ml, PP
------	-------------------------

Product number	72.689.004
----------------	------------

Product description	Reaction tube, working volume: 2 ml, material: PP, transparent, cap: yellow, retaining loop cap, cap attached, 500 piece(s)/bag
---------------------	---

Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

Color of product	transparent
------------------	-------------

Closure material	Low Density Polyethylene (LD-PE)
------------------	----------------------------------

Color of cap	yellow
--------------	--------

Size

Diameter	10.8 mm
----------	---------

Length of product	39 mm
-------------------	-------

Product characteristics

Graduation	no
------------	----

Cap	Retaining loop cap
-----	--------------------

Centrifugation max (RCF)	20000 x g
--------------------------	-----------

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	5000
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	150000
Height of case	340 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0539 cbm
Weight of case	7.48 kg
EAN inner box	4038917116459
EAN case	4038917015141

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT