



## Product specification

### Reaction tube, 2 ml, PP



## Product information

Product number	72.691.004
----------------	------------

## Product descriptions

Name	Reaction tube, 2 ml, PP
Product number	72.691.004
Product description	Reaction tube, working volume: 2 ml, material: PP, yellow, cap: yellow, retaining loop cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag

## Material & colours

Color of product	yellow
Closure material	Polypropylene (PP)
Color of cap	yellow

## Size

Diameter	10.8 mm
Length of product	40 mm

## Product characteristics

Graduation	yes
Cap	Retaining loop cap
Closure type	Hinged lid
Centrifugation max (RCF)	20000 x g

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	120000
Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	6.05 kg
EAN inner box	4038917312882
EAN case	4038917312875

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT