



Product specification

Screw cap micro tube, 1.5 ml



Product information

Product number	72.692.100
----------------	------------

Product descriptions

Name	Screw cap micro tube, 1.5 ml
Product number	72.692.100
Product description	Screw cap micro tube, working volume: 1.5 ml, conical base, with knurling, transparent, cap: natural, cap attached, 500 piece(s)/bag

Material & colours

Color of product	transparent
O-Ring material	Ethylene propylene diene rubber (EPDM)
Closure material	Polypropylene (PP)
Color of cap	natural

Size

Diameter	10.8 mm
Length of product	44 mm

Product characteristics

Base shape	Conical base
Graduation	no
Knurls	yes

Cap	Retaining loop cap
Closure type	Screw cap
Centrifugation max (RCF)	20000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	150000
Height of case	340 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0539 cbm
Weight of case	8.12 kg
EAN inner box	4038917116855
EAN case	4038917015363

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT