



## Product specification

### Transport case I 50, internal case for block rack



## Product information

Product number	95.1720.011
----------------	-------------

## Product descriptions

Name	Transport case I 50, internal case for block rack
Product number	95.1720.011
Product description	Transport case I 50, internal case for block rack, (LxWxH): 270 x 170 x 215 mm

## Material & colours

Product material	Acrylnitrile butadiene styrene (ABS)
------------------	--------------------------------------

## Size

Width of product	170 mm
Length of product	270 mm
Height of product	215 mm

## Product characteristics

Design	internal case for block rack
Shape of inner tube	square

## Purity & certification

Satisfies the requirement	TRBA 100
Product category	no medical device   IVD

Batched	no
---------	----

## Packaging

Piece(s) / pallet	108
Height of case	195 mm
Width of case	223 mm
Depth of case	308 mm
Volume of case	0.0134 cbm
Weight of case	1.34 kg
EAN inner box	4038917159951
EAN case	4038917159944

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT