



Product specification

Battery bag, (LxW): 205 x 115 mm, PL



Product information

Product number	95.1725
----------------	---------

Product descriptions

Name	Battery bag, (LxW): 205 x 115 mm, PL
Product number	95.1725
Product description	Battery bag, without batteries, (LxW): 205 x 115 mm, material: PL, for standard cooling batteries, 12 piece(s)/case

Material & colours

Product material	polyester (PL)
------------------	----------------

Size

Width of product	115 mm
Length of product	205 mm

Product characteristics

Design	without batteries
--------	-------------------

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order quantity	12
Piece(s) / inner box	12
Piece(s) / case	12
Piece(s) / pallet	1980
Height of case	120 mm
Width of case	228 mm
Depth of case	338 mm
Volume of case	0.0092 cbm
Weight of case	5.5 kg
EAN inner box	4038917087568
EAN case	4038917068376

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT