

Product specification

Screw cap micro tubes, 2 ml, sterile



Product information

Order number	72.694.011
--------------	------------

Product descriptions

Name	Screw cap micro tubes, 2 ml, sterile
Product number	72.694.011
Product description	Screw cap micro tubes, working volume: 2 ml, skirted conical base, with knurling, transparent, cap: yellow, cap assembled, sterile

Product characteristics

Base shape	skirted conical base
Graduation	no
Knurls	yes
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	44 mm

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Color of cap	yellow

Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation

Batched

yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	96000
Height of case	140 mm
Width of case	278 mm
Depth of case	348 mm
Volume of case	0.0135 cbm
Weight of case	1.86 kg
EAN inner box	4038917117210
EAN case	4038917067744

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

