



## Product specification

### Screw cap micro tube, 2 ml



## Product information

Product number	72.694.100
----------------	------------

## Product descriptions

Name	Screw cap micro tube, 2 ml
Product number	72.694.100
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, cap: natural, cap attached, 500 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
O-Ring material	Ethylene propylene diene rubber (EPDM)
Closure material	Polypropylene (PP)
Color of cap	natural

## Size

Diameter	10.8 mm
Length of product	44 mm

## Product characteristics

Base shape	skirted conical base
------------	----------------------

Graduation	no
Knurls	yes
Cap	Retaining loop cap
Closure type	Screw cap
Centrifugation max (RCF)	20000 x g

## Purity & certification

Satisfies the requirement	IATA
Product category	diagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order quantity	5000
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	120000
Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	9.06 kg
EAN inner box	4038917117593
EAN case	4038917015677

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT