



Product specification

Screw cap micro tube, 2 ml



Product information

| | |
|----------------|------------|
| Product number | 72.694.100 |
|----------------|------------|

Product descriptions

| | |
|---------------------|--|
| Name | Screw cap micro tube, 2 ml |
| Product number | 72.694.100 |
| Product description | Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, cap: natural, cap attached, 500 piece(s)/bag |

Material & colours

| | |
|------------------|--|
| Color of product | transparent |
| O-Ring material | Ethylene propylene diene rubber (EPDM) |
| Closure material | Polypropylene (PP) |
| Color of cap | natural |

Size

| | |
|-------------------|---------|
| Diameter | 10.8 mm |
| Length of product | 44 mm |

Product characteristics

| | |
|------------|----------------------|
| Base shape | skirted conical base |
| Graduation | no |

| | |
|--------------------------|--------------------|
| Knurls | yes |
| Cap | Retaining loop cap |
| Closure type | Screw cap |
| Centrifugation max (RCF) | 20000 x g |

Purity & certification

| | |
|---------------------------|--------------------------------------|
| Satisfies the requirement | IATA |
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Batched | yes |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 500 |
| Piece(s) / case | 5000 |
| Piece(s) / pallet | 120000 |
| Height of case | 415 mm |
| Width of case | 398 mm |
| Depth of case | 398 mm |
| Volume of case | 0.0657 cbm |
| Weight of case | 9.06 kg |
| EAN inner box | 4038917117593 |
| EAN case | 4038917015677 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT