



Product specification

Screw cap micro tube, 2 ml, sterile



Product information

Product number	72.694.106
----------------	------------

Product descriptions

Name	Screw cap micro tube, 2 ml, sterile
Product number	72.694.106
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, cap: natural, cap attached assembled, with printed writing space, with graduations, sterile, 100 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
O-Ring material	Ethylene propylene diene rubber (EPDM)
Closure material	Polypropylene (PP)
Color of cap	natural

Size

Diameter	10.8 mm
Length of product	44 mm

Product characteristics

Base shape	skirted conical base
------------	----------------------

Label/Print	With print
Color of print/label	white
Graduation	yes
Knurls	yes
Cap	Retaining loop cap
Centrifugation max (RCF)	20000 x g

Purity & certification

Satisfies the requirement	IATA
Product category	diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	96000
Height of case	140 mm
Width of case	278 mm
Depth of case	348 mm
Volume of case	0.0135 cbm
Weight of case	1.98 kg
EAN inner box	4038917117647
EAN case	4038917015691

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT