



# Product specification

## Screw cap micro tube, 2 ml



### Product information

Product number	72.694.107
----------------	------------

### Product descriptions

Name	Screw cap micro tube, 2 ml
Product number	72.694.107
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, cap: natural, cap attached, with printed writing space, with graduations, 500 piece(s)/bag

### Material & colours

Color of product	transparent
O-Ring material	Ethylene propylene diene rubber (EPDM)
Closure material	Polypropylene (PP)
Color of cap	natural

### Size

Diameter	10.8 mm
Length of product	44 mm

### Product characteristics

Base shape	skirted conical base
Label/Print	With print

Color of print/label	white
Graduation	yes
Knurls	yes
Cap	Retaining loop cap
Centrifugation max (RCF)	20000 x g

## Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	120000
Height of case	415 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0657 cbm
Weight of case	9.07 kg
EAN inner box	4038917117654
EAN case	4038917015707

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT