



Product specification

Transport case I 100, internal case with attachment



Product information

Product number	95.1720
----------------	---------

Product descriptions

Name	Transport case I 100, internal case with attachment
Product number	95.1720
Product description	Transport case I 100, internal case with attachment, (LxWxH): 360 x 310 x 275 mm

Size

Width of product	310 mm
Length of product	360 mm
Height of product	275 mm

Product characteristics

Design	internal case with attachment
Shape of inner tube	square

Purity & certification

Satisfies the requirement	TRBA 100
Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / pallet	30
Height of inner box	50 mm
Width of inner box	90 mm
Depth of inner box	337 mm
Height of case	325 mm
Width of case	358 mm
Depth of case	398 mm
Volume of case	0.0463 cbm
Weight of case	2.05 kg
EAN inner box	4038917087551
EAN case	4038917076210

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT