



Product specification

Transport case B 13, flat design for block rack



Product information

| | |
|----------------|-------------|
| Product number | 95.1717.011 |
|----------------|-------------|

Product descriptions

| | |
|---------------------|--|
| Name | Transport case B 13, flat design for block rack |
| Product number | 95.1717.011 |
| Product description | Transport case B 13, flat design for block rack, (LxWxH): 395 x 395 x 300 mm |

Material & colours

| | |
|------------------|-----------------------------------|
| Product material | High Density Polyethylene (HD-PE) |
|------------------|-----------------------------------|

Size

| | |
|-------------------|--------|
| Width of product | 395 mm |
| Length of product | 395 mm |
| Height of product | 300 mm |

Product characteristics

| | |
|--------|----------------------------|
| Design | flat design for block rack |
|--------|----------------------------|

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | P 650, ADR |
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------|---------------|
| Piece(s) / pallet | 30 |
| Height of case | 340 mm |
| Width of case | 398 mm |
| Depth of case | 398 mm |
| Volume of case | 0.0539 cbm |
| Weight of case | 3.3 kg |
| EAN inner box | 4038917133197 |
| EAN case | 4038917069250 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT