



## Product specification

### SafeSeal reaction tube, 2 ml, PP, PCR Performance Tested



## Product information

Product number	72.695.400
----------------	------------

## Product descriptions

Name	SafeSeal reaction tube, 2 ml, PP, PCR Performance Tested
Product number	72.695.400
Product description	SafeSeal reaction tube, working volume: 2 ml, material: PP, transparent, cap: transparent, SafeSeal cap, cap attached, with injected graduation and writing space, PCR Performance Tested, 250 piece(s)/minigrip bag

## Material & colours

Color of product	transparent
Closure material	Polypropylene (PP)
Color of cap	transparent

## Size

Diameter	10.8 mm
Length of product	40 mm

## Product characteristics

Graduation	yes
Cap	SafeSeal cap

Centrifugation max (RCF)	25000 x g
--------------------------	-----------

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Piece(s) / inner box	1000
Piece(s) / case	2000
Piece(s) / pallet	100000
Height of inner box	161 mm
Width of inner box	221 mm
Depth of inner box	301 mm
Height of case	180 mm
Width of case	313 mm
Depth of case	458 mm
Volume of case	0.0258 cbm
Weight of case	3.15 kg
EAN inner box	4038917223188
EAN case	4038917223171

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT