



Product specification

Reaction tube, 0.5 ml, PP



Product information

Product number	72.699.003
----------------	------------

Product descriptions

Name	Reaction tube, 0.5 ml, PP
Product number	72.699.003
Product description	Reaction tube, working volume: 0.5 ml, material: PP, blue, cap: blue, standard cap, cap attached, 500 piece(s)/bag

Material & colours

Color of product	blue
Closure material	Polypropylene (PP)
Color of cap	blue

Size

Diameter	7.8 mm
Length of product	30 mm

Product characteristics

Graduation	no
Cap	Standard cap
Closure type	Hinged lid
Centrifugation max (RCF)	20000 x g

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	320000
Height of case	215 mm
Width of case	298 mm
Depth of case	398 mm
Volume of case	0.0255 cbm
Weight of case	2.87 kg
EAN inner box	4038917118156
EAN case	4038917015905

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT