



Product specification

Closure clip, (LxW): 220 x 12 mm, PA



Product information

Product number	95.1714
----------------	---------

Product descriptions

Name	Closure clip, (LxW): 220 x 12 mm, PA
Product number	95.1714
Product description	Closure clip, (LxW): 220 x 12 mm, material: PA, for cross-bottom bag 95.1713, 25 piece(s)/case

Material & colours

Product material	Polyamide (PA)
------------------	----------------

Size

Width of product	12 mm
Length of product	220 mm

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order quantity	25
Piece(s) / inner box	25
Piece(s) / case	25

Piece(s) / pallet	8100
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	1.4 kg
EAN inner box	4038917087520
EAN case	4038917065948

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT