



## Product specification

### Reaction tube, 0.5 ml, PP



## Product information

Product number	72.699.004
----------------	------------

## Product descriptions

Name	Reaction tube, 0.5 ml, PP
Product number	72.699.004
Product description	Reaction tube, working volume: 0.5 ml, material: PP, green, cap: green, standard cap, cap attached, 500 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	green
Closure material	Polypropylene (PP)
Color of cap	green

## Size

Diameter	7.8 mm
Length of product	30 mm

## Product characteristics

Graduation	no
Cap	Standard cap
Closure type	Hinged lid
Centrifugation max (RCF)	20000 x g

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order quantity	5000
Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	320000
Height of case	215 mm
Width of case	298 mm
Depth of case	398 mm
Volume of case	0.0255 cbm
Weight of case	2.87 kg
EAN inner box	4038917118163
EAN case	4038917015912

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT