



Product specification

Screw cap micro tube, 1.5 ml, Biosphere® plus



Product information

| | |
|----------------|------------|
| Product number | 72.703.217 |
|----------------|------------|

Product descriptions

| | |
|---------------------|--|
| Name | Screw cap micro tube, 1.5 ml, Biosphere® plus |
| Product number | 72.703.217 |
| Product description | Screw cap micro tube, working volume: 1.5 ml, skirted conical base, with knurling, transparent, cap: natural, cap attached, with printed writing space, with graduations, Biosphere® plus, 25 piece(s)/bag |

Material & colours

| | |
|------------------|--|
| Color of product | transparent |
| O-Ring material | Ethylene propylene diene rubber (EPDM) |
| Closure material | Polypropylene (PP) |
| Color of cap | natural |

Size

| | |
|-------------------|---------|
| Diameter | 10.8 mm |
| Length of product | 44 mm |

Product characteristics

| | |
|------------|----------------------|
| Base shape | skirted conical base |
|------------|----------------------|

| | |
|--------------------------|--------------------|
| Label/Print | With print |
| Color of print/label | white |
| Graduation | yes |
| Knurls | yes |
| Cap | Retaining loop cap |
| Centrifugation max (RCF) | 20000 x g |

Purity & certification

| | |
|---------------------------|--------------------------------------|
| Satisfies the requirement | IATA |
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Purity standard | Biosphere plus |
| Sterilisation | Gas sterilisation |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 250 |
| Piece(s) / case | 500 |
| Piece(s) / pallet | 49500 |
| Height of inner box | 141 mm |
| Width of inner box | 191 mm |
| Depth of inner box | 211 mm |
| Height of case | 160 mm |
| Width of case | 223 mm |
| Depth of case | 398 mm |
| Volume of case | 0.0142 cbm |
| Weight of case | 1.45 kg |
| EAN inner box | 4038917223836 |
| EAN case | 4038917223829 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

