



## Product specification

**Screw cap micro tube, 1.5 ml, PCR Performance Tested,  
Low protein-binding**



### Product information

Product number	72.703.600
----------------	------------

### Product descriptions

Name	Screw cap micro tube, 1.5 ml, PCR Performance Tested, Low protein-binding
Product number	72.703.600
Product description	Screw cap micro tube, working volume: 1.5 ml, skirted conical base, with knurling, transparent, cap: natural, cap assembled, with printed writing space, with graduations, PCR Performance Tested, Low protein-binding, 100 piece(s)/minigrip bag

### Material & colours

Color of product	transparent
O-Ring material	Ethylene propylene diene rubber (EPDM)
Closure material	Polypropylene (PP)
Color of cap	natural

### Size

Diameter	10.8 mm
Length of product	44 mm

Max. volume	1.5 ml
-------------	--------

## Product characteristics

Base shape	skirted conical base
Label/Print	With print
Color of print/label	white
Graduation	yes
Knurls	yes
Cap	standard screw cap
Centrifugation max (RCF)	20000 x g

## Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	1000
Piece(s) / pallet	99000
Height of inner box	141 mm
Width of inner box	191 mm
Depth of inner box	211 mm
Height of case	160 mm
Width of case	223 mm
Depth of case	398 mm
Volume of case	0.0142 cbm
Weight of case	2.17 kg
EAN inner box	4038917393089
EAN case	4038917393072

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support  
and is valid without signature.



# SARSTEDT