



Product specification

SafeSeal reaction tube, 0.5 ml, PP



Product information

Product number	72.704
----------------	--------

Product descriptions

Name	SafeSeal reaction tube, 0.5 ml, PP
Product number	72.704
Product description	SafeSeal reaction tube, working volume: 0.5 ml, material: PP, transparent, cap: transparent, SafeSeal cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	Polypropylene (PP)
Color of cap	transparent

Size

Diameter	7.8 mm
Length of product	30 mm

Product characteristics

Graduation	yes
Cap	SafeSeal cap

Centrifugation max (RCF)	30000 x g
--------------------------	-----------

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	2000
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	156000
Height of inner box	97 mm
Width of inner box	151 mm
Depth of inner box	189 mm
Height of case	135 mm
Width of case	358 mm
Depth of case	388 mm
Volume of case	0.0188 cbm
Weight of case	1.61 kg
EAN inner box	4038917162180
EAN case	4038917162173

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT