




## Product specification

### SafeSeal reaction tube, 0.5 ml, PP



## Product information

|                |            |
|----------------|------------|
| Product number | 72.704.002 |
|----------------|------------|

## Product descriptions

|                     |   |
|---------------------|---|
| Name                | SafeSeal reaction tube, 0.5 ml, PP  |
| Product number      | 72.704.002  |
| Product description | SafeSeal reaction tube, working volume: 0.5 ml, material: PP, red, cap: red, SafeSeal cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag |

## Material & colours

|                  |                    |
|------------------|--------------------|
| Product material | Polypropylene (PP) |
| Color of product | red                |
| Closure material | Polypropylene (PP) |
| Color of cap     | red                |

## Size

|                   |        |
|-------------------|--------|
| Diameter          | 7.8 mm |
| Length of product | 30 mm  |

## Product characteristics

|                          |              |
|--------------------------|--------------|
| Graduation               | yes          |
| Cap                      | SafeSeal cap |
| Centrifugation max (RCF) | 30000 x g    |

## Purity & certification

|                  |                       |
|------------------|-----------------------|
| Product category | diagnostic device, CE |
| CE certified     | CE                    |
| Batched          | yes                   |

## Packaging

|                        |               |
|------------------------|---------------|
| Minimum order quantity | 2000          |
| Piece(s) / inner box   | 500           |
| Piece(s) / case        | 2000          |
| Piece(s) / pallet      | 156000        |
| Height of inner box    | 97 mm         |
| Width of inner box     | 151 mm        |
| Depth of inner box     | 189 mm        |
| Height of case         | 135 mm        |
| Width of case          | 358 mm        |
| Depth of case          | 388 mm        |
| Volume of case         | 0.0188 cbm    |
| Weight of case         | 1.61 kg       |
| EAN inner box          | 4038917227551 |
| EAN case               | 4038917227544 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT