



Product specification

Cross-bottom bag, (LxW): 270 x 270 mm, PE



Product information

Product number	95.1713
----------------	---------

Product descriptions

Name	Cross-bottom bag, (LxW): 270 x 270 mm, PE
Product number	95.1713
Product description	Cross-bottom bag, (LxW): 270 x 270 mm, material: PE, for transport case B 17, 25 piece(s)/case

Size

Width of product	270 mm
Length of product	270 mm

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / inner box	25
Piece(s) / case	25
Piece(s) / pallet	1800
Height of inner box	28 mm

Width of inner box	90 mm
Depth of inner box	165 mm
Height of case	125 mm
Width of case	303 mm
Depth of case	333 mm
Volume of case	0.0126 cbm
Weight of case	1.7 kg
EAN inner box	4038917087513
EAN case	4038917065917

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT