



## Product specification

### SafeSeal reaction tube, 0.5 ml, PP



## Product information

Product number	72.704.006
----------------	------------

## Product descriptions

Name	SafeSeal reaction tube, 0.5 ml, PP
Product number	72.704.006
Product description	SafeSeal reaction tube, working volume: 0.5 ml, material: PP, blue, cap: blue, SafeSeal cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag

## Material & colours

Color of product	blue
Closure material	Polypropylene (PP)
Color of cap	blue

## Size

Diameter	7.8 mm
Length of product	30 mm

## Product characteristics

Graduation	yes
Cap	SafeSeal cap
Centrifugation max (RCF)	30000 x g

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	156000
Height of inner box	97 mm
Width of inner box	151 mm
Depth of inner box	189 mm
Height of case	135 mm
Width of case	358 mm
Depth of case	388 mm
Volume of case	0.0188 cbm
Weight of case	1.61 kg
EAN inner box	4038917227643
EAN case	4038917227636

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT