



## Product specification

**SafeSeal reaction tube, 0.5 ml, PP, PCR Performance**

**Tested, Low protein-binding**



## Product information

Product number	72.704.600
----------------	------------

## Product descriptions

Name	SafeSeal reaction tube, 0.5 ml, PP, PCR Performance Tested, Low protein-binding
Product number	72.704.600
Product description	SafeSeal reaction tube, working volume: 0.5 ml, material: PP, transparent, cap: transparent, SafeSeal cap, cap attached, with injected graduation and writing space, PCR Performance Tested, Low protein-binding, 50 piece(s)/minigrip bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	Polypropylene (PP)
Color of cap	transparent

## Size

Diameter	7.8 mm
Length of product	30 mm

## Product characteristics

Graduation	yes
Cap	SafeSeal cap
Centrifugation max (RCF)	30000 x g

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Minimum order quantity	1200
Piece(s) / inner box	300
Piece(s) / case	1200
Piece(s) / pallet	93600
Height of case	135 mm
Width of case	358 mm
Depth of case	388 mm
Volume of case	0.0188 cbm
Weight of case	1.2804 kg
EAN inner box	4038917265508
EAN case	4038917265492

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT