



Product specification

SafeSeal reaction tube, 1.5 ml, PP



Product information

Product number	72.706.001
----------------	------------

Product descriptions

Name	SafeSeal reaction tube, 1.5 ml, PP
Product number	72.706.001
Product description	SafeSeal reaction tube, working volume: 1.5 ml, material: PP, brown, cap: brown, SafeSeal cap, cap attached, with injected graduation and writing space, 250 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	brown
Closure material	Polypropylene (PP)
Color of cap	brown

Size

Diameter	10.8 mm
Length of product	39 mm

Product characteristics

Graduation	yes
Cap	SafeSeal cap
Centrifugation max (RCF)	30000 x g

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	2000
Piece(s) / inner box	1000
Piece(s) / case	2000
Piece(s) / pallet	100000
Height of inner box	161 mm
Width of inner box	221 mm
Depth of inner box	301 mm
Height of case	180 mm
Width of case	313 mm
Depth of case	458 mm
Volume of case	0.0258 cbm
Weight of case	2.96 kg
EAN inner box	4038917118224
EAN case	4038917078245

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT