



Product specification

SafeSeal reaction tube, 1.5 ml, PP, Biosphere® plus



Product information

Product number	72.706.200
----------------	------------

Product descriptions

Name	SafeSeal reaction tube, 1.5 ml, PP, Biosphere® plus
Product number	72.706.200
Product description	SafeSeal reaction tube, working volume: 1.5 ml, material: PP, transparent, cap: transparent, SafeSeal cap, cap attached, with injected graduation and writing space, Biosphere® plus, 50 piece(s)/bag

Material & colours

Color of product	transparent
Closure material	Polypropylene (PP)
Color of cap	transparent

Size

Diameter	10.8 mm
Length of product	39 mm

Product characteristics

Graduation	yes
Cap	SafeSeal cap
Centrifugation max (RCF)	30000 x g

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	250
Piece(s) / case	500
Piece(s) / pallet	72000
Height of inner box	126 mm
Width of inner box	151 mm
Depth of inner box	201 mm
Height of case	145 mm
Width of case	213 mm
Depth of case	323 mm
Volume of case	0.01 cbm
Weight of case	0.9 kg
EAN inner box	4038917118279
EAN case	4038917077491

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT