

# Product specification

## Screw cap micro tubes, 0.5 ml



### Product information

Order number	72.730.003
--------------	------------

### Product descriptions

Name	Screw cap micro tubes, 0.5 ml
Product number	72.730.003
Product description	Screw cap micro tubes, working volume: 0.5 ml, skirted conical base, with knurling, transparent, without cap, 500 piece(s)/bag

### Product characteristics

Base shape	skirted conical base
Graduation	no
Knurls	yes
Centrifugation max (RCF)	20000 x g

### Size

Diameter	10.8 mm
Length of product	44 mm

### Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

### Purity & certification

Satisfies the requirement	IATA
Accordance with MDD or IVDD	In-vitro diagnostic, CE
CE certified	CE
Lot or serial number	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	180000
Height of case	300 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0475 cbm
Weight of case	6.52 kg
EAN inner box	4038917119252
EAN case	4038917016520

**This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.**

**This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.**

This document was prepared by EDP support and is valid without signature.

